



REACTHEALTH

## New Luna G3 CPAP and Auto-CPAP



Luna G3 CPAP  
SKU: LG3500

Luna G3 Auto-CPAP  
SKU: LG3600

React Health's all new Luna G3 CPAP and Auto-CPAP devices come standard with all the features that make sense for you and your patients. Our standard integrated cellular modem makes capturing patient compliance a breeze.

Our standard patient-forward features include integrated heated tubing, pre-heat humidifier function and Auto Ramp offering added comfort to patient therapy as well as a bedside clock and resupply reminders for the convenience of the patient. You will find these features and more, in a smaller footprint, with an easy-to-read full color screen and intuitive interface.

# Device Specifications

Device Dimensions (L x W x H)	10.4" x 5.7" x 4.5"
Weight	3.75 lbs.
Mode of Operation	CPAP (LG3500), CPAP, Auto-CPAP (LG3600)
SD Card	Standard
Pressure Range	4 to 20cmH2O (in 0.5 cm increments)
Ramp Time	0 to 60 minutes
Sound	<26 dBA when device is at 10cmH2O
Water Capacity	350 ml (to maximum fill line)

## Replacement Parts

LG 3500

LG 3600

LG 34510

LG 3HT

UN 100

LG31020\*

LG380088X\*

LG3B

LPG3

LPAC

PRTG301

Luna G3 CPAP Device - HT

Luna G3 APAP Device - HT

Luna G3 Replacement Water Chamber

Luna G3 Integrated Heated Tubing

3B Ultra Noir CPAP Standard Tubing (Non-Heated)

Luna G3 Disposable Air Filter, 2 per bag, 3 per box (White)

Luna G3 Reusable Filter - Box of 12 Filters (Black)

Luna G3 Carry Bag

Luna G3 Power Supply Adapter

Luna G2/G3 Power Cord

Luna G3 Repl Filter Cover- Black



REACTHEALTH

[www.reacthealth.com](http://www.reacthealth.com)

[info@reacthealth.com](mailto:info@reacthealth.com)